



DV2000 Diagnostic Imaging Software Supported Image Sources

Image Capture Source	Capture Type	Windows 7 (32bit)	Windows 7 (64 bit)	Windows 8.1 (64 bit)	Windows 10 (64 bit)
TWAIN	Still Image	No ¹	No ¹	No	No
Leutron	Video	Yes	Yes	Yes	Yes
Piccolo ³	Video	No	No	No	No
Direct Show	Video	Yes	Yes	Yes	Yes
WIA	Still Image	No ¹	No ¹	No	No
Canon EOS RC ³	Still Image	No	No	No	No
Canon EOS ED	Still Image	Compatibility depends on driver support which varies between camera models. Please consult Canon website www.canon.com			
Nikon D100 ³	Still Image	No	No	No	No
AVT Marlin	Video	Yes	Yes	Yes	Yes
Artray	Video	No ¹	No ¹	No ¹	No ¹
Elite / Alkeria	Video	Yes	No	No	No
File Watcher	Still Image	Yes	Yes	Yes	Yes
File Import	Still Image	Yes	Yes	Yes	Yes
Point Grey	Video	Yes ²	No ¹	Yes	Yes
Keeler	Video	Yes	Yes	Yes	Yes
CSO	Video	Yes	Yes	Yes	Yes
Inami	Video	Refer to Point Grey image source compatibility			
Sun Kingdom	Video	Refer to Canon EOS ED image source compatibility			
Flashbus ³	Video	No	No	No	No
E300 Video Converter	Video	Yes	Yes	Yes	Yes

¹ Not tested. This image source may be used on this operating system at customer discretion but no support will be offered if it functions incorrectly.

² Some Point grey cameras have been shown to perform poorly on computers with low system memory, in these situations reducing the camera frame rate to 30 fps or less will provide better performance.

³ This image source is no longer supported and is only available for those who insist on running Studio on unsupported operating system. These installations will not be supported by Medmont and availability of this image source is not guaranteed in future releases.